

# The Center for Health Promotion and Preventive Medicine (USACHPPM) Gas Chromatographic Procedures for Analysis of Water and Soil for Energetics and Related Compounds

Mike Hable and Curtis Oliver

# Topics of Discussion

- USACHPPM Method
- Other Methods Available for the Detection and Measurement of Explosives
- Advantages and Disadvantages of Methodologies

# USACHPPM Extraction Method


## ■ Water Extraction

- 100 mL sample
- 1.0 mL IAA
- Rotate for 30 Min


## ■ Soil Extraction

- 2-5g sample, 20 mL H<sub>2</sub>O
- 5.0 mL IAA
- Shake for 3 Hr

# Rotary Agitators



# Flatbed Shaker




14 3 2003

# Method Reporting and Detection Limits

<u>Compound</u>	<u>Water - µg/L</u>		<u>Soil - µg/g*</u>	
	<u>MDL</u>	<u>MRL</u>	<u>MDL</u>	<u>MRL</u>
2,6-Dinitrotoluene	0.001	0.01	0.004	0.01
2,4-Dinitrotoluene	0.003	0.02	0.006	0.02
2,4,6-Trinitrotoluene	0.004	0.03	0.007	0.01
RDX	0.016	0.1	0.006	0.01
HMX	0.58	3	0.023	0.05
2-Nitrotoluene	0.016	0.09	0.006	0.02
3-Nitrotoluene	0.044	0.09	0.005	0.02
4-Nitrotoluene	0.035	0.09	0.006	0.02
Nitrobenzene	0.004	0.03	0.005	0.02
1,3-Dinitrobenzene	0.011	0.09	0.005	0.02
1,3,5-Trinitrobenzene	0.004	0.03	0.006	0.02
4-Amino-2,6-Dinitrotoluene	0.023	0.1	0.026	0.05
2-Amino-4,6-Dinitrotoluene	0.026	0.1	0.009	0.02
Tetryl	0.023	0.5	0.012	0.02
Nitroglycerin	0.015	0.09	0.025	0.05

\* Soil data based on a 5.0-gram sample.

# Chromatogram



# EPA Method 8330/8095

- Liquid Chromatography (8330)

- Water
  - Soil

- Gas Chromatography (8095)

- Water
  - Soil

# Method Reporting and Detection Limits EPA 8330

Soil - µg/g\*

<u>Compound</u>	<u>MDL</u>	<u>MRL</u>
2,6-Dinitrotoluene	0.1	0.25
2,4-Dinitrotoluene	0.1	0.25
2,4,6-Trinitrotoluene	0.1	0.25
RDX	0.4	1
HMX	0.5	2.2
2-Nitrotoluene	0.1	0.25
3-Nitrotoluene	0.1	0.25
4-Nitrotoluene	0.1	0.25
Nitrobenzene	0.1	0.25
1,3-Dinitrobenzene	0.1	0.25
1,3,5-Trinitrobenzene	0.1	0.25
4-Amino-2,6-Dinitrotoluene	0.1	0.25
2-Amino-4,6-Dinitrotoluene	0.1	0.25
Tetryl	0.1	0.65
Nitroglycerin	—	—


- Soil data based on a 2.0-gram sample.

# Chromatogram

Page 1 of 1

Software Version : 6.2.0.0.0:S27  
Sample Name : 2 ppm 8330  
Instrument Name : MAXIMUS  
Rack/Vial : 1/1  
Sample Amount : 1.000000  
Cycle : 1  
Date : 8/7/02 2:30:10 PM  
Data Acquisition Time : 8/7/02 1:43:44 PM  
Channel : A  
Operator : bishoppw  
Dilution Factor : 1.000000

Result File :  
Sequence File : \wg6099\MUSCP\RWBISHOP\LCdata\test.seq



## Explosives in Soil

Peak #	Time [min]	Area [ $\mu$ V·s]	Height [ $\mu$ V]	Component Name	ug/mL
6	4.616	189155.14	9665.50	HMX	1.3714
7	7.249	252717.86	8582.68	RDX	1.3512
9	10.100	518041.28	20458.30	1,3,5TNB	1.3012
12	13.110	720240.15	24550.58	1,3DNB	1.1504
13	13.731	386876.86	12044.66	Tetryl	1.5899
14	15.609	517497.29	15153.44	NB	1.3858
15	16.752	498602.27	16407.29	2,4,6TNT	1.1949
16	17.638	223607.36	6374.85	4A26DNT	0.7961
17	18.646	416149.34	9743.19	2A46DNT	0.9707
18	20.361	356410.76	8970.48	2,6DNT	1.2504
19	22.070	711026.71	16356.39	2,4DNT	1.2514
23	26.330	296136.40	6690.56	2NT	1.0972
25	29.556	291132.90	5875.41	4NT	1.2360
26	31.435	370058.88	7140.17	3NT	1.2566

5747653.20 168013.49

# Advantages

## ■ USACHPPM Method

- Simple, reproducible
- Mini-Extraction
- High sample throughput
- Low detection limits
- Applicable to deployment samples

# Disadvantages

- Purity of isoamyl acetate
- Emulsions

# USACHPPM Method

## ■ Additional Applications

- Vegetation Samples
- Blood
- Destruction of Munitions (PMACWA)
- Range Residue
- Air Sampling

## ■ Future Applications

- Biological tissue