



Toxicogenomics :

Understanding the molecular consequences
of low-level toxicant exposure

U.S. Army SBCCOM

Edgewood Chemical Biological Center

Principal Project Scientist:

Jennifer Weeks Sekowski, Ph.D.

Molecular Engineering Team Leader:

James J. Valdes, Ph.D.

Collaborators: Mary Anne Orehek, Stephanie Mahmoudi, Jennifer Bucher, Robert Mioduszewski, Christopher Whalley, Bernard Benton, Stan Hulet, Sandra Thomson, Maryanne Vahey and Martin Nau



Why is Gene Expression Important?

Edgewood Chemical Biological Center




Plants



Animals



In vitro cell culture




Toxicants,
drugs, and other factors
can influence
the copying of the gene
message

DNA – the genetic code

Messenger RNA (mRNA)-
the working copy of a gene

**Too little, too much, or
incorrectly coded
gene product**
can disrupt normal cellular
processes leading to injury and
disease in the organism

Gene product – what the cell
makes from the mRNA







How DNA microarrays (gene chips) work

Edgewood Chemical Biological Center

Sample mRNA

Experimental

Control



Mix and hybridize cDNA

DNA Microarray



Spot (gene) 1 2 3

DNA Microarray



Spot (gene)


1

2

3

Fluorescence
Detection

Laser Excitation




Spot (gene) 1 2 3



Bioinformatics to analyze the raw data,

Edgewood Chemical Biological Center




GeneSpring™ software screenshot of gene expression levels measured on an Affymetrix Rat Neurobiology Gene Chip®



... find genes with similar expression profiles,

Edgewood Chemical Biological Center

Expression level



Time

GeneSpring™
software screen
shot of gene
expression
levels measured
on an Affymetrix
Yeast Gene
Chip®



... identify the functional class of a gene,

Edgewood Chemical Biological Center



GeneSpring™
software
screenshot of
gene functional
classes.

Gene
expression was
measured on an
Affymetrix Yeast
Gene Chip®




...and tell us what's known about the gene

Edgewood Chemical Biological Center

Link to web based databases :

- 1) Genbank (NCBI)- Gene sequence data
- 2) Swiss-Prot (ExPASy) – Protein data
- 3) PubMed (NCBI), ToxNet (NLM)- Literature information



Hear Res 1994 Aug;78(2):235-42

Identification of a glutamate/aspartate transporter in the rat cochlea.


Li HS, Niedzielski AS, Beisel KW, Hiel H, Wenthold RJ, Morley BJ

The neurotransmitter at the synapses between hair cells and spiral ganglion cells in the cochlea is probably L-glutamate or a similar excitatory amino acid. Glutamate uptake by nerve terminals and glial cells is an important component of neurotransmission at glutamatergic synapses of the central nervous system, for providing a reservoir of transmitter or transmitter precursors and the termination of the released glutamate. Hair cell synapses are not surrounded by glial cells, therefore, the uptake mechanism for glutamate in the cochlea may be unique. cDNA was synthesized from total RNA isolated separately from the rat organ of Corti, spiral ganglia, and lateral wall tissues. The expression of a glutamate/aspartate transporter (GLAST) was detected by DNA amplification with the polymerase chain reaction. The other two members of glutamate transporters in this family were not detected by this method....




Scheme of a Toxicogenomics Study

Edgewood Chemical Biological Center



mRNA extracted

fluorescently-labeled cDNA
hybridized to gene array



read out of expression levels

array analysis bioinformatics
(e.g. GeneSpring™ software & web
based databases)

Validation by
correlation with
other
toxicological
analyses

- Pathology
- Histology
- Physiology
- Enzymology
- Behavior

Could the gene product play a role in
injury or the initiation or promotion of a
disease ?

Changes at the gene level can precede
disease by years

Confirm results by RT-PCR

Link altered
genes to
**known
functions**

- Metabolism
- Cell cycle regulation
- Neurotransmitter
- Inflammatory response



Applications of Gene Array Analysis

Edgewood Chemical Biological Center

Toxicological and Medical Investigation

- Finding biomarkers for earlier diagnosis of disease
- Finding early warning (“sentinel”) biomarkers for toxicity
- Identifying molecular “signatures” of drug or chemical exposure
- Identifying molecular targets for therapeutic intervention
- Potential to ID genes that predispose individuals to chemical/drug sensitivities
- Advantage in high throughput screening of drugs and chemicals
- Helps in the reduction of animal use




Applications of Gene Array Analysis

Edgewood Chemical Biological Center


Biological Agent Detection

Unknown sample



Extract mRNA from sample,

Convert to fluorescently-labeled cDNA,
Hybridize cDNA to Bioagent Detection
Microarray




Fluorescence
Reader



Bioinformatics
translates raw
data

Detectors



Warfighter Model

Clinician Model

- 1) B.anthraxis, strain X
- 2) B. anthracis, strain Y**
- 3) B. anthracis, strain Z
- 4) B. abortus
- 5) B. melitensis**
- 6) B. suis
- 7) E.coli (EPEC)
- 8) E.coli (EHEC)**




Our Gene Array Facility can...

Edgewood Chemical Biological Center

Affymetrix Gene Chip® Custom DNA Microarray

Genes arrayed here



Genes arrayed here

1. Prepare the RNA sample, convert to labeled cDNA, hybridize
2. Carry out all steps necessary to measure the fluorescent signal from the DNA microarray.
3. Determine relative gene expression levels in control and experiment
4. Determine genes with altered expression
5. Analyze altered genes of interest via literature search and laboratory experimentation