Early and Rapid Diagnosis of Exposure to Biothreat Agents Using Host Gene Expression Responses

Rina Das and Marti Jett

Dept of Mol Pathology, WRAIR

Pathogen specialists / collaborators....

Rina, Pas, Roger, Neill, George Ludwig, R. Hammamieh, C. Mendis, Marti, Jett

WRAIR

- David Hoover (Brucella)
- Luther Lindler (Plague)
- Robert Putnak (Dengue)

USAMRIID (Anthrax, VEE, Botulinum toxin, & Dengue viruses)

- George Ludwig
- Erik Henchal
- Melody Jensen & Len Smith
- Nirenjan Kanessa-Thassan

Univ MD-Biotechnology Center

Dr. Maria Salvato (influenza, Arena viruses, Lassa)

■ WRAIR-Mol Path (SE toxins, LPS)

Rina Das Roger Neill
Chanaka Mendis, Chris Cummings,
Preveen Ramamoorthy, S. Bi
Apsara Dhokalia,
Sachin Mani

Rasha Hammamieh-Bioinformatics Biostaticians from Mitre Corp BioSpice project

Georgetown- Chemistry

Dr. David Yang (Cholera toxin)

Atabak Royaee

Walter Reed Army Institute of Research

HOST GENE EXPRESSION RESPONSES TO BIOLOGICAL THREAT AGENTS:

- -Rapidand Early Diagnosis of Pathogen Exposure
- -Design of Stage-specific Therapeutic Approaches

Traditional approach

- Direct pathogen identification methods are limited due to:
- -dependence on concentration/ time
- -tissue sequestration of pathogen
- Mutants (natural or deliberate)
 Chimeras, etc, can elude detection
 by structural-based probes or other
 direct methods

Host response approach

- Peripheral blood lymphoid cells circulate throughout (minutes):
 - -cells encounter pathogens
 - -they respond and create a record of the induced "insults".
 - Draw blood sample, 1-2ml, automatically isolate RNA, gene profile indicates pathogen exposure & stage of illness




Create a library of host gene expression responses to biological threat agents

- using peripheral blood mononuclear cells (PBMC) as the target tissue, isolate RNA, perform microarray analysis on these samples.
- Broadly compare major biological threat agents for multiple exposure time periods in vitro to PBMC from healthy donors
- Confirm in vivo using monkeys exposed to the threat agent and take PBMC at 3-4 time periods (3 NHP per time period).



TECHNOLOGY INVOLVED:

- Use library of ~40,000 cDNA gene fragments
 - Spotted onto 3 slides in duplicate




The "test" is compared to a reference sample in order to provide maximum continuity over time.

Control and test are separately compared to the reference.

Reference is a mixture of tumor samples that permits normalization over years of study.



Pathogen induced Gene Profiles




These inflammatory mediators may not differentiate among pathogenic agents, but may still be useful markers to gauge illness progression

Principal Component Analysis

A method for analyzing total gene pattern to illustrate the differences observed for various toxins.


PCA Mapping of s2 (81.3%)



To identify differentiating gene profiles, algorithms were applied


- Statistical power evaluated externally
 - 2-4 experiments for each of 3-5 exposure times per pathogenic agent
- Bonferroni multiple test correction applied to strengthen the reliability of pattern identification

SPECIFIC GENE PATTERNS FOR EACH PATHOGEN





Confirmation of cDNA microarray data




Selection of a few genes that differentiate Brucella exposure from the 8 other pathogens studied

Gene responses as a function of time post-exposure

Non-human primate gene patterns
24, 48, 72 h post exposure to B. anthracis
at 8 LD₅₀


(mild flu-like illness began ~72 h)

Gene Clustering Analysis in blood samples of NHP exposed to *B anthracis*



OTSG Casualty Projections LTC Deborah Schnelle

Projected survival rates vs time treatment is initiated




Analysis platform

- The 40,000 gene chip is a splendid research tool but not a viable option in the field
- We envision selecting the "quiet" genes that we have characterized as a basis for diagnosis (~100-1000)
- An electrochemical, microfluidics chip has some of the features we seek-Nanogen
- Metrigenix 4D array platform is another promising technology




Advantages of Nanogen's Microelectronic Array

- Electronic movement and deposition
- Electronic concentration
- Electronic hybridization



The Nanogen gene expression platform could be a very attractive tool in translational research as in tracking biomarkers in clinical samples. It offers the advantage of multiplexing many targets.

Metrigenix Technology



The patented Flow-Thru Chip technology (FTC) serves as the technology base for MetriGenix. The FTC is a state-of-the-art 4-dimensional microarray technology platform that can be used to provide high-throughput, high-content assays to measure gene expression activity. The advantages of developing Flow-Thru Chips include:

- •improved assay responsivity due to the increased surface area
- •reduced assay times due to enhanced mass-transport within the channels
- •more uniform spot morphology due to the wetting properties of the microporous chip
- •smaller sample and reagent volume requirements due to the reduction in the reaction volume




Host Gene Expression Profiles Induced by Biological Threat Agents:

Provides rapid/early diagnosis plus tailored therapeutic approaches

Rapidly compare the gene expression profiles (1000's of genes) with a library of known biological threat and other pathogenic agents

Determine gene pattern to learn: Extent of Individual Exposure, Individual Susceptibility, Course of Impending Illness and Undetectable Potentiating Contaminants

Draw ~2 ml of blood from an exposed individual, use automated system to isolate RNA.



Gene groups

A B C D E F1 F2 G H I J

Anthrax
BoTox A
Brucella
Dengue
Plague
SEB
VEE

Rapid analysis of exposure. continue monitoring the "negative" profiles

Initiate treatments based on course of impending illness using pattern for agents with similar gene expression profiles.

Standard Dipstick toxin detection method does/does not determine the agent present. Chaos and panic are eminent



Gene Array Technology: Summary

Uses for Gene Array Technology

- Attack by an unknown agent
- Outbreak of an unknown illness

Specific Diagnosis

- Causative Agent
- Degree of Progression of Illness

Treatment

- Mechanism of pathologic action of toxin or infection
- Suggests specific treatments

■ The Future

- Electrochemical detection on microchips
- Microfluidics and 4D arrays

